

Improving the reliability of in vitro toxicity and ecotoxicity measurements with nanomaterials

Elijah J. Petersen, Shannon Hanna, Monique Johnson, Gregory Cooksey, Bryan Nelson, John Elliott, Lee Yu


**National Institute of Standards and Technology
(NIST)**

Presented 3/10/2015

Overview

- 1. Introduction/Background**
- 2. Use of cause-and-effect analysis to design a robust nanoecotoxicity assay with *C. elegans***
- 3. Single-particle ICP-MS method to quantify the size distribution of gold nanoparticle uptake by *C. elegans***

NIST Role in Nano-EHS





National Nanotechnology Initiative 2011 Environmental Health and Safety Research Strategy

NIST & NanoToxicology

Three Major Areas

- Standards
- Tools for Nanomaterial Characterization
- Methods/Data to inform health and environmental risk models
- Robust methods to enable reliable toxicity testing
- Metrology for evaluating transport and transformations of nanomaterials in biological and environmental systems



NIST Standard Reference Materials

Gold nanoparticles (10, 30, and 60 nm)

Single-wall carbon nanotube (raw soot) and dispersed into three length populations

Titanium dioxide nanoparticles (made from Degussa P25)

2 nm silicon nanoparticles

Silver nanoparticle (75 nm PVP coated)

Silver nanoparticles (10 nm in preparation)

Multiwall carbon nanotube (in preparation)

Can be useful for interlaboratory comparisons, instrument validation and calibration, and positive and negative controls for nanotoxicity studies

Critical for establishing comparability of nano-related measurements.



Documentary Standards



NIST participates in standards organizations that provide validated documentary standards on a range of topics

- Nanoparticle characterization using a range of instruments for all nanoparticles (DLS, TEM, etc.) through the NIST/NCL protocols
- Sonication protocols that provide reproducible, traceable NP sonication between instruments and laboratories
- MTS assay for cell toxicity from nanomaterials
- Guidance document for aquatic toxicity testing of nanomaterials

Overview

1. **Introduction/Background**
2. **Use of cause-and-effect analysis to design a robust nanoecotoxicity assay with *C. elegans***
3. **Single-particle ICP-MS method to quantify the size distribution of gold nanoparticle uptake by *C. elegans***


Why *C. elegans*?

- Abundant in soils and sediments
- > 70% of *C. elegans* proteome has human homologues
- 12 out of 17 signal transduction pathways conserved between *C. elegans* and humans
- Complete cell lineage map, large knockout libraries, established genetic methods
- “easy” to culture
- Small size & short lifespan
- Higher-throughput *in vivo* assays
- Transparent




ISO Method 10872

Preparation




Assay






- Uses positive control benzylcetyldimethylammonium chloride (BAC C16 – $EC_{50} = 15.1 \text{ mg l}^{-1}$)
- Only test specification is growth inhibition of 20-80% at EC_{50} value

Automated *C. elegans* Imaging



C. elegans Image Stitching



This is an open access article published under an ACS AuthorChoice [License](#), which permits copying and redistribution of the article or any adaptations for non-commercial purposes.



Perspective

pubs.acs.org/crt


Chemical Research in Toxicology

Use of Cause-and-Effect Analysis to Design a High-Quality Nanocytotoxicology Assay


Matthias Rösslein,^{†,⊥} John T. Elliott,^{*,§,⊥} Marc Salit,^{||} Elijah J. Petersen,[§] Cordula Hirsch,[†] Harald F. Krug,[‡] and Peter Wick[†]

Poster: Toward achieving harmonization in a nano-cytotoxicity assay measurement through an interlaboratory comparison study

Cause & Effect Analysis of *C. elegans* Assay



EC₅₀ of BAC C16 Dose-Response Curve



Growth inhibition of 15mg/L BAC C16







Plate shaking




Influence of media on BAC toxicity



Influence of live/dead feed on bioassay




Bioassay feed dependence



Tests were conducted with 15 mg/L BAC-C16

Impact of worm culture on bioassay










*This experiment was performed with the plates shaking

Plate design to test for NP artifacts


Plates 1-3

Plate 4




-  Negative control
-  Positive chemical control
-  Test sample
-  ENM control #1
-  ENM control #2

Poster: Identification and avoidance of potential artifacts and misinterpretations in nanomaterial ecotoxicity measurements



Polystyrene Nanoparticle Toxicity Assay




EC50s for growth

- ◆ 12/12/2014 100.7
- 12/18/2014 125
- △ 12/29/2014 81.7
- × 1/2/2015 74.9
- * 1/9/2015 119.6

Overview

- 1. Introduction/Background**
- 2. Usage of cause-and-effect analysis to design a robust nanoecotoxicity assay with *C. elegans***
- 3. Single-particle ICP-MS method to quantify the size distribution of gold nanoparticle uptake by *C. elegans***

Single Particle ICP-MS (spICP-MS)




- Operates ICP-MS in time resolved analysis mode: a particle event lasts 0.4 ms, typical dwell time 3 - 10 ms
- Particle size distribution and number concentration

* *J. Anal. At. Spectrom.*, 2012, 27, 1143

Base Digestion of *C. elegans* for spICP-MS analysis


- 7 % TMAH (w/v) used for dissolution of worm tissue
- Dilution to bring TMAH concentration to <0.1%
- sp-ICP-MS measurement to test NPs accumulation on the organisms

spICP-MS Measurements of control *C. elegans*




Single particle ICP-MS signals for *C. elegans* in control exposure.
Mean Au intensity: 0.18 counts. Date of exposure: June 24, 2014

spICP-MS Measurements of 30 nm AuNPs after digestion procedure (no worms)




Single particle ICP-MS signals for 30 nm AuNPs (NIST) following treatment with 7% TMAH. **Control digestion, no worms.**

Size distribution of 30 nm AuNPs (NIST) extracted following TMAH treatment. **Control digestion, no worms.**




spICP-MS Measurements of *C. elegans* exposed to 30 nm AuNPs




Single particle ICP-MS signals for *C. elegans* spiked with 30 nm AuNPs (NIST). 799 particle events in 60 s. **Diluted by a factor of 1500.**

Size distribution of NPs extracted from *C. elegans* exposed to 30 nm AuNPs (NIST).




spICP-MS Measurements of 60 nm AuNPs after digestion procedure (no worms)




Single particle ICP-MS signals for 60 nm AuNPs (NIST) following treatment with 7% TMAH. Control digestion, no worms.

Size distribution of 60 nm AuNPs (NIST) extracted following TMAH treatment. Control digestion, no worms.




spICP-MS Measurements of *C. elegans* exposed to 60 nm AuNPs







Single particle ICP-MS signals for *C. elegans* spiked with 60 nm AuNPs (NIST). 412 particle events in 60 s. Diluted by a factor of 2500.

Size distribution of NPs extracted from *C. elegans* exposed to 60 nm AuNPs (NIST).







Sucrose density gradient centrifugation to separate *C. elegans* and AuNPs







Spin in
centrifuge





% Recovery of Au conc: $95.6\% \pm 3.7\%$



% Recovery: 95.6% ± 3.7%

Conclusions

1. Cause and effect analysis was used to identify major sources of uncertainty in a *C. elegans* standard method
2. Plate shaking and bacteria concentration were shown to have the strongest impact on assay results
3. Growth inhibition variability was similar for BAC-C16 and PSNPs
4. A method was developed to analyze AuNP size distributions in *C. elegans*

Cause & Effect Analysis: A new approach for developing robust bionano assays



St. Gallen, Switzerland

June 18 & 19, 2015

[http://www.empa.ch/plugin/template/empa/22/
155222/---/l=2](http://www.empa.ch/plugin/template/empa/22/155222/---/l=2)